

## Investigation of retention mechanisms in HILIC chromatography: Important considerations for robust method development



Anders Fridström  
Sigma- Aldrich GMBH

© 2010 Sigma-Aldrich Co. All rights reserved.

**SIGMA-ALDRICH**

### Abstract

#### Investigation of Retention Mechanisms in HILIC Chromatography: Important Considerations for Robust Method Development

Anders Fridström<sup>1</sup>, David S. Bell<sup>2</sup>

<sup>1</sup>Sigma-Aldrich, Buchs SG/Switzerland, <sup>2</sup>Supelco/Sigma-Aldrich, Bellefonte/USA  
Anders.fridstrom@sial.com

Hydrophilic interaction liquid chromatography (HILIC), especially in conjunction with mass spectrometry (MS), has become a powerful tool for the analysis of a wide variety of challenging analytes. Applications of the technique have increased dramatically over the past decade, especially for the analysis of polar analytes where reversed-phase chromatography suffers. HILIC conditions employ a high percentage of acetonitrile which enables facilitated solvent evaporation in LC/MS sources and thus often an increase in analyte response when compared to more aqueous based systems. The increased retention of polar analytes afforded by HILIC provides improved selectivity and decreases the impact of endogenous species, often leading to improved qualitative and quantitative analyses [1]. Although HILIC has proven useful, it has also been thwarted with complications including difficulties in method development and method robustness.

In this presentation, studies investigating the underlying retention mechanisms dominant in HILIC chromatography are presented and discussed. Along with reversed-partitioning HILIC is well known to exhibit, ion-exchange and the interplay of the dominant mechanisms are unveiled and used to develop a model of overall retention and selectivity. Interactions that operate using different stationary phase chemistries and conditions are presented. The impact of analyte polarity and charge as well as the variations caused by high percentages of organic on these physicochemical parameters are highlighted. Throughout the discussion, examples of use and misuse of HILIC are employed to illustrate these important concepts to build a solid fundamental foundation for efficient and effective use of this powerful technique.

© 2010 Sigma-Aldrich Co. All rights reserved.

**SIGMA-ALDRICH**

## Agenda



Introduction

Factors affecting the HILIC Separation


Modelling Ionic interactions on polar stationary phases in HILIC

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Introduction

HILIC Hydrophilic Interaction Liquid Chromatography



© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Factors effecting the HILIC systems

Column

Mobile Phase

- pH
- Buffer Concentration

Analytes

- pKa
- Log P<sub>OW</sub> or Log D<sub>OW</sub>

(Temperature)

Sample

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Columns

**Ascentis Express**  
Extreme Performance on Any LC System


Ascentis Express Fused Core Particle Columns

HILIC (Si)

OH5 (Pentalol) Branched hydroxylated alkane

New!!


F5 (PFP)



© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## pH Effect of Acetonitrile on pH of Ammonium Acetate [4]



Measurements were taken at 25°C.

Triangle: 90.0% ACN,  
Square: 75% ACN,  
Diamond: 50% ACN,  
Circle: 32.5% ACN

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## pKa

Analyte  $pK_a$  values have also been shown to be impacted by the presence of organic modifiers


Figure 2 shows the results of an NMR experiment conducted that explored the chemical shift of a proton near the ionizable group for amitriptyline in 90% acetonitrile. From data such as this, effective  $pK_a$  values can be established for a variety of compounds.

Table 1 shows the results for several basic pharmaceutical compounds. The data indicates that the effective  $pK_a$  value for a basic analyte in 90% acetonitrile is approximately 1  $pK_a$  unit less than the aqueous-based value

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Determination of $pK_a$ Values using $^1H$ NMR [4]



© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Table 1: Determination of $pK_a$ Values using NMR [4]

Amitriptyline			Analyte	Literature $pK_a$	$pK_a_w^s$	Correlation ( $R^2$ )
% Acetonitrile	$pK_a$	Correlation ( $R^2$ )				
25	9.32	0.9997	Amitriptyline	9.4	8.34	0.9923
50	9.02	0.9996	Nortriptyline	9.7	8.92	0.9920
75	8.88	0.9956	Diphenhydramine	9.0	8.33	0.9978
90	8.34	0.9923	Verapamil	8.9	7.98	0.9976
			Alprenolol	9.7	8.73	0.9855

- $pK_a$  values for bases decrease with increasing acetonitrile
- At 90% each analyte exhibited a  $pK_a$  value about 1 full pH unit less than the literature  $pK_a$  value

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Investigation of Retention Mechanisms on Different HILIC Phases

interaction differences for three different HILIC stationary phases: PFPP (F5), bare silica (HILIC) and a new pentalol phase (OH5)

Using ephedrine as a probe molecule, retention as a function of buffer concentration was collected and interpreted.




Related compounds and dominant interactions prevalent using each phase is studied

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Selected Probes

ACD/Labs, PhysChemProp, v. 12

Structure	pKa(MB)	LogD(8.0)	LogP	MW	name
	9.38	-0.37	1.08	165.23	pseudoephedrine
	9.38	-0.37	1.08	165.23	ephedrine
	9.37	-1.35	0.13	167.21	synephrine

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Experimental

Conditions:

Instrument: #9, Waters 2690/Micromass ZQ single quadrupole interfaced via ESI operating in pos. ion mode

Columns:

- Ascentis Express Pentanol (OH5), 10 cm x 3.0 mm
- Ascentis Express HILIC, 10 cm c 3.0 mm,
- Ascentis Express F5, 10 cm x 3.0 mm,

Mobile Phase A: **10 mM ammonium acetate** (pH unadjusted) in 10:90 water:acetonitrile

Mobile Phase B: 10:90 water:acetonitrile

**Mixtures of A and B were run at 0%B, 20%B, 40%B, 60%B and 80%B corresponding to 10 mM, 8 mM, 2 mM, 4 mM and 2 mM buffer concentrations, respectively**

Flow rate: 0.4 mL/min

Temperature: ambient

Detection: MS, ESI, pos ion mode, scan m/z 125 – 300

Injection volume: 2 uL

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Calculation for ion exchange impact in HILIC

Samples were injected in triplicate at each buffer concentration using each of the phases. A sample of ephedrine only was also injected under each condition to discriminate from pseudoephedrine in the mix.

$$\text{Log } k = -\text{log}[C^+]_m + \text{log}\beta_{\text{IEX}}$$

$[C^+]_m$  concentration of the competing ion in the mobile phase and

$\beta_{\text{IEX}}$  constant for a given system

- phase ratio,
- ion-exchange capacity of the stationary phase ion-exchange equilibrium constant.

Log k Retention


A plot of log k vs log  $[C^+]_m$  will thus yield a slope of -1 when ion-exchange is solely responsible for retention,

The plot would yield a slope of 0 where ion-exchange is not present.

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®


## Response of Ephedrine Retention on Buffer Concentration on Three HILIC Phases



© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## HILIC OH5 column at 10 and 2 mM ammonium acetate




© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®



### HILIC (bare silica) column at 4 and 2 mM ammonium acetate

Syneprine still last  
More Rt change  
Both partition and IEX




© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®


### F5 column at 8 and 2 mM ammonium acetate

Syneprine early!  
Rt change with Buffer  
Little or no partition,  
IEX present



© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®



## Conclusions


Columns described in this presentation

- Pentalol (OH5)– dominated by partition
- Bare silica (HILIC)– both partition and IEX
- Pentafluorophenylpropyl – mainly IEX

In order to develop robust and reliable methods using HILIC chromatography the following factors should be considered

- pH changes in high % organics
- pKa of analytes in high % organics
- Mobile phase modifiers concentrations

## Method Development and Optimization Chart



21

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## References


1. Hydrophilic Interaction Liquid Chromatography (HILIC) and Advanced Applications, Wang Perry G., He Weixuan, CRC Press, Taylor & Francis Group.
2. Needham, S.R., Bell, D., J. Chromatogr., A. 2000, 869, 159-170.
3. McCalley, D. V., J. Chromatogr., A. 2010, 1217, 3408-3417.
4. Dinh, N. P., Jonsson T., Irgum K., J. Chromatogr., A. 2010, 1217, 3408-3417.
5. W. Naidong, Journal of Chromatography B 796 (2003) 209.
6. D.S. Bell, Jones, A. Daniel, Journal of Chromatography A 1073 (2005) 99.
7. D.S. Bell, Brandes, Hillel K., in 30th International Symposium and Exhibit on High Performance Liquid Phase Separations and Related Techniques, San Francisco, California USA, 2006.
8. D.S. Bell, Solute Attributes and Molecular Interactions Contributing to Retention on a Fluorinated High-Performance Liquid Chromatography Stationary Phase, Thesis, The Pennsylvania State University, 2005
9. ACD PhysChem, v. 12, Advanced Chemistry Development, Toronto, ON Canada

© 2010 Sigma-Aldrich Co. All rights reserved.

SIGMA-ALDRICH®

## Acknowledgements

- Dave Bell
- Hugh Cramer
- Craig Aurand
- Wayne Way
- Gaurang Parmar



23

© 2010 Sigma-Aldrich Co. All rights reserved.

**SIGMA-ALDRICH**

# Thanks!

© 2010 Sigma-Aldrich Co. All rights reserved.

**SIGMA-ALDRICH**